



226 Creekstone Ridge
Woodstock, GA 30188

MODEL NCU5400 SHARPS KIOSK - PRODUCT SPECIFICATIONS

Model NCU5400 is a self-standing heavy-duty, locking unit in which the waste cannot be accessed by the public once the sharps are deposited.

- Features:**
- sharp kiosk dimension – 54 “ x 23” x 24”
 - unit accommodates up to an 18 gallon US sharps container
 - pull out drawer to facilitate removing the container
 - optional white interior floor for better visibility
 - heavy 14 gauge steel construction
 - powder coated steel
 - drop box opening to fit 2 quarts containers
 - disclaimer inside door which can be customized to district and state laws
 - customer color available with large orders
 - manufactured in US

Security Features

- dual locking system for high safety
- drop box opening prevents public from having access to sharps
- anti-pry lower door
- unit able to be anchored to the ground
- funnel system has been designed to direct used sharps directly into
- selected container

Labelling Specifications

- signage dimension – 18” x 24 “
- (2) contrasting colors with reflective lettering
- legible from a 10 foot distance
- (1) inch lettering for header (4) line of text
- half inch lettering for remainder of text
- wording - custom made as per client requirement

Optional Sharps Disposal US Government approved Mail Back System

- specially designed mail back sharps container
- postage prepaid return shipping box with bar code tracking
- a protective 4–mil bag liner
- instructions in both English and Spanish
- a multi–part simplified manifest tracking form

**WHETHER YOU USE A MAIL BACK PROGRAM OR AN IN HOUSE
PICK UP SERVICE – OUR SYSTEM IS DESIGNED FOR BOTH**

